







# Mindset24 Software Reference Manual



**NOLAN**  
NOLAN COMPUTER SYSTEMS, L.L.C.

# ***MINDSET24 SOFTWARE REFERENCE MANUAL***

Software Version 1.0

for use with the

**MINDSET Model Mindset24**

**NOLAN COMPUTER SYSTEMS, L.L.C. (NCS)**

**Mindset Software Users Manual**

**Software Version 1.0 July 2003**

Copyright ' 2001 - 2003 Nolan Computer Systems, L.L.C. All Rights Reserved.

The information in this manual is furnished for informational use only and subject to change without notice. Nolan Computer Systems assumes no responsibility or liability for any errors or inaccuracies that may appear within.

Printed in the USA.  
July 2003

## IMPORTANT - PLEASE READ CAREFULLY

This End-User License Agreement ( License ) is a legal agreement between you and NOLAN COMPUTER SYSTEMS, L.L.C. ( NCS ), the manufacturer. This License covers the MINDSET software, the associated media, any printed materials and any on-line or electronic documentation. Any software provided by third parties is covered by a separate warranty. You may use this NCS SOFTWARE only in connection with the use of the MINDSET.

**HIGH RISK ACTIVITIES:** MINDSET is not fault-tolerant and is not designed, manufactured or intended for use or resale as on-line control equipment in hazardous environments requiring fail-safe performance, such as in the operation of nuclear facilities, aircraft navigation or communication systems, air traffic control, direct life support machines or weapons systems, in which the failure of MINDSET could lead directly to death, personal injury, or severe physical or environmental damage ( High Risk Activities ). NCS and its distributors specifically disclaim any expressed or implied warranty of fitness for High Risk Activities. MINDSET is a research instrument and is not intended for use in clinical medical applications and is not Food and Drug Administration (FDA) approved. NCS and its distributors specifically disclaim any expressed or implied warranty of fitness for medical or therapeutic use.

**NCS SOFTWARE LICENSE:** The NCS SOFTWARE is protected by copyright laws and international copyright treaties, as well as other intellectual property laws and treaties. NCS SOFTWARE is licensed, not sold.

**GRANT OF LICENSE:** This License grants you the following rights:

- You may install and use one copy of the NCS SOFTWARE on your computer.
- **Back-up Copy:** you may, and should, create a back-up copy of the NCS SOFTWARE included with MINDSET. You may use the back-up copy solely for archival purposes.

**DESCRIPTION OF OTHER RIGHTS AND LIMITATIONS:**

- **Limitation on Reverse Engineering, Decompilation and Disassembly.** You may not modify, reverse engineer, decompile, or disassemble the NCS SOFTWARE in whole or in part.
- **Single MINDSET.** The NCS SOFTWARE is licensed with the MINDSET as a single integrated product. The NCS SOFTWARE may only be used with MINDSET.
- **Rental.** You may not rent or lease the NCS SOFTWARE.
- **Software Transfer.** You may permanently transfer all of your rights under this License only as part of a sale or transfer of MINDSET, provided you retain no copies, transfer all of the NCS SOFTWARE (including all copies, component parts, the media and printed materials, all versions and any upgrades of the NCS SOFTWARE and this License), and the recipient agrees to the terms of this License.
- **Termination.** Without prejudice to any other rights, NCS may terminate this License if you fail to comply with the terms and conditions of this License. In such an event, you must destroy all copies of the NCS SOFTWARE and all of its component parts.

**COPYRIGHT:** All title and copyrights in and to the NCS SOFTWARE (including but not limited to any images, photographs, animation, video, audio, music, text and applets incorporated into the NCS SOFTWARE) and any copies of the NCS SOFTWARE are owned by NCS. All rights not specifically granted under this License are reserved by NCS.

**U.S. GOVERNMENT RESTRICTED RIGHTS:** The software product and documentation are provided with RESTRICTED RIGHTS. Use, duplication, or disclosure by the United States Government is subject to restrictions as set forth in subparagraph (c)(i) and (2) of the Commercial Computer Software-Restricted Rights at 48 CFR 52.227-19, as applicable.

## SOFTWARE WARRANTY

**LIMITED MEDIA WARRANTY:** NCS warrants that for a period of ninety (90) days from the date of its delivery to you the media on which the NCS SOFTWARE is furnished will be free from defects in materials and workmanship under normal use. This limited warranty extends only to you as the original licensee. NCS's entire liability and your exclusive remedy will be replacement of the media meeting NCS's limited warranty which is returned to NCS with proof of purchase in the form of a bill of sale as evidence that the media is within the warranty period. NCS will have no responsibility to replace a distribution disc/disk damaged by accident, abuse or misapplication.

Any implied warranties on the media, including the implied warranties of merchantability and fitness for a particular purpose, are limited in duration to ninety (90) days from the date of delivery. Some states do not allow limitations on how long an implied warranty lasts, so these limitations may not apply to you. This warranty gives you specific legal rights and you may also have other rights which vary from state-to-state and country-to-country.

**EXCLUSION OF WARRANTY ON NCS SOFTWARE:** You expressly acknowledge and agree that use of the NCS SOFTWARE is at your sole risk. The NCS SOFTWARE is provided as is and without warranty of any kind and NCS and NCS's licensors (hereinafter, NCS and NCS's licensors shall be collectively referred to as NCS) expressly disclaim all warranties, expressed or implied, including, but not limited to, the implied warranties of merchantability and fitness for a particular purpose. NCS does not warrant that the functions contained in the NCS SOFTWARE will meet your requirements or that the operation of the NCS SOFTWARE will be corrected. Furthermore, NCS does not warrant or make any representations regarding the use or the results of the use of the NCS SOFTWARE in terms of its correctness, accuracy, reliability or otherwise. No oral or written information or advice given by NCS or an NCS authorized representative shall create a warranty or in any way increase the scope of this warranty. Should the NCS SOFTWARE prove defective, you (and not NCS or an NCS authorized representative) assume the entire cost of all necessary servicing, repair or correction. Some jurisdictions do not allow the exclusion of implied warranties, so the above exclusion may not apply to you.

**LIMITATION OF LIABILITY:** Hereinafter, NCS and NCS's licensors shall be collectively referred to as NCS. NCS shall not be liable for any incidental or consequential damages for breach of any expressed or implied warranty, breach of contract, negligence, strict liability or any other legal theory related to this product. Such damages include, but are not limited to, loss of profits, loss of revenue, loss of data, loss of use of the product or any associated equipment, down time and purchaser's time, even if NCS has been advised of the possibility of such damages. In any case, NCS's entire liability under any provision of this agreement shall be limited to the amount actually paid by you allocable to the NCS software. Some states do not allow the exclusion or limitation of consequential or incidental damages, so the above exclusion or limitation may not apply to you.

Some states do not allow the exclusion or limitation of incidental or consequential damages, or allow limitations on how long an implied warranty lasts, so the above limitations or exclusions may not apply to you. This warranty gives you specific legal rights and you may have other rights which vary from state-to-state and country-to-country.

## **RESPONSIBILITY OF THE USER**

This product and its components can work properly only if it is implemented and maintained as provided by the instructions in this manual and the hardware manual for the Mindset24. In case of any failure, the product should not be used until completely repaired by an authorized distributor. The owner of this product is responsible for any failure in the operation caused by ill use or inappropriate maintenance, or due to repairs made by persons not authorized by NCS.

Users are responsible for the implementation and use of this product.

### **NOTE**

Copyrighted (©), trademarked (™) and registered (®) products are mentioned in this document. Rather than indicate that the product/software is copyrighted, trademarked or registered each time the product/software is mentioned, the authors of this document hereby state that there is no intent to infringe on said copyright, trademark or registration.

# Table of Contents

1	INTRODUCTION.....	1
1.1	Assumptions and Conventions of This Manual .....	1
2	SOFTWARE INSTALLATION.....	3
2.1	Minimum Requirements .....	3
2.2	Preparation and Installation.....	3
2.3	Starting MindMeld .....	3
2.4	System Check-out.....	4
3	MINDMELD CONCEPTS .....	5
3.1	Basic MindMeld Concepts.....	5
3.1.1	Streaming Data .....	5
3.1.2	Array Data.....	5
4	THE MINDMELD24 LIVE DATA CAPTURE APPLICATION .....	6
4.1	MindMeld24 Live Data Capture Main Window .....	6
4.1.1	MindMeld24 Live Data Capture Main Window Controls .....	6
4.2	MindMeld24 Live Data Capture Tools Overview .....	10
4.3	Tools In Detail .....	12
4.3.1	Mindset MS24 Tool.....	12
4.3.2	Storage Tool .....	15
4.3.3	Stripchart Tool.....	17
4.3.4	Fast Fourier Transform (FFT) Tool .....	21
4.3.5	Averager Tool .....	25
4.3.6	Histogram Tool.....	26
4.3.7	Neuromapper Tool .....	29
4.3.8	Compressed Spectral Array (CSA) Tool .....	31
4.3.9	Hemispheric Tool .....	35
5	CREATE/EDIT/RUN A MINDMELD LIVE DATA CAPTURE SCHEME.....	37
5.1	Using the Create/Edit Scheme Window .....	37
5.2	Running a Scheme.....	39
6	THE MINDMELD ANALYSIS FUNCTIONS APPLICATION .....	40
6.1	MindMeld Analysis Functions Main Window Controls .....	40
6.2	MindMeld Analysis Functions: Create/Edit/Replay Scheme.....	46
6.2.1	MindMeld24 Analysis Functions Mindset 1 Tool Properties .....	47
6.2.2	MindMeld24 Analysis Functions Statistics Window.....	49
6.2.3	MindMeld24 Analysis Functions Storage Tool .....	51
6.2.4	MindMeld24 Analysis Functions Stripchart Tool .....	51
6.2.5	MindMeld24 Analysis Functions Remaining Tools.....	51
	Appendix A: Acronyms and Glossary .....	A-1
	Appendix B: Using the Mindset Calibrator .....	B-1
	Appendix C: Listing of Pre-configured Schemes .....	C-1

# MINDSET SOFTWARE USERS MANUAL

## 1 INTRODUCTION

Welcome to the Mindset24 application software known as MindMeld24. The Mindset hardware and MindMeld24 software combine to provide you with powerful electroencephalography (EEG) data capture and analysis techniques. These capabilities were formerly available only to hospitals, universities and other deep pocketed institutions due to their complexity and expense.

The MindMeld24 software is divided into two separate sub-applications:

- *MindMeld24 Live Data Capture (LDC)* is used for data acquisition and real-time data analysis. MindMeld24 captures and analyzes live EEG data. Various tools within MindMeld24 are implemented by a *scheme* methodology used for selecting and connecting the tools together. One tool allows you to save data to your hard drive for off-line analysis. The sample pre-configured schemes illustrate proper tool configurations. You can design and save new schemes to focus on your particular area of interest.
- *MindMeld24 Analysis Functions (AF)* is used for off-line analysis of a patient's session, also known as *replay*. After the EEG data has been captured and saved to disk with MindMeld24 LDC, the tools in MindMeld24 AF can be applied to the data to provide an overall interpretation of the patient's session. One very useful option in MindMeld24 AF allows you to remove artifacts and other invalid data prior to sending the data through the analysis tools, thereby improving diagnostic accuracy.

### 1.1 ASSUMPTIONS AND CONVENTIONS OF THIS MANUAL

This manual assumes you are familiar with the operation of your Personal Computer (PC) and your Windows operating system (O/S). You should know how to insert a Compact Disc — Read-only Memory (CD-ROM) disc and start an application program, use your mouse to select objects on the screen and to implement the basic capabilities of a graphical user interface (GUI). If you have questions about these basic operations, refer to your PC's hardware manuals, Windows manuals or on-line Help.

This manual uses the term *System* to mean the combination of the MindMeld24 applications, an electro-cap, and the Mindset24 hardware appropriately configured to acquire and analyze EEGs. The Mindset Calibrator hardware is an optional hardware/software instrument for measuring the frequency alignment of your Mindset. Calibrator details are in Appendix B.

*This typeface is used to draw your attention to information that is of particular importance.*

This symbol → is used to mean the GUI process of selecting from the main menu thru the required option; for example, *Scheme → Edit Scheme*.

This manual does not attempt to cover all the dialog boxes and information/error messages that appear as you use MindMeld24. Follow the on-screen prompts to acknowledge the message/dialog or to complete the operation.

EEG information shown in the various MindMeld24 windows in this document is for example purposes only.

**NOTE:** There are some comments and graphics in this document that allude to capabilities that are not included in this release. These capabilities will be present in a later software release.

## 2 SOFTWARE INSTALLATION

This section instructs you on how to install the MindMeld24 software which is distributed on CD-ROM. Floppy disks are available upon request.

### 2.1 MINIMUM REQUIREMENTS

Before beginning, ensure that your PC meets the following minimum requirements:

- a Pentium class or equivalent central processing unit (CPU). For optimum real-time performance, a *200 Mega-Hertz (MHz) Pentium (or faster) is recommended.*
- Windows 95 or later O/S, or Windows NT 4.0 or later O/S.
- at least 32 megabytes of Random Access Memory (RAM).
- at least 20 megabytes of available hard drive space.
- at least a video graphics adapter (VGA) class video adapter. An accelerated graphics adapter greatly enhances real-time neuromapping display performance.
- a Windows compliant Small Computer Systems Interface (SCSI) adapter and cable.
- a printer installed with a Windows printer driver.

### 2.2 PREPARATION AND INSTALLATION

Ensure that there are no other application programs running on your PC. If there are, first save any work in these applications and then quit them.

Insert the CD-ROM disc into your CD-ROM drive, the MindMeld24 installer program should start automatically. If not, open the icon representing your CD-ROM drive and double-click on the install file. Then follow the on-screen instructions.

Once installation is completed, a directory named MindMeld24 is created on your hard drive containing the following items:

- the MindMeld24 Live Data Capture (LDC) and MindMeld24 Analysis Functions (AF) applications
- Sample Scheme subdirectory containing some sample pre-configured schemes (.em3 files)
- Sample Data subdirectory containing some sample EEG files.

All of the required software is now on your PC and ready to run.

### 2.3 STARTING MINDMELD24

Start MindMeld24 by:

- navigating to the MindMeld24 folder and double-clicking your required MindMeld24 application:
  - ♦ MindMeld24.exe for MindMeld24 LDC
  - or*
  - ♦ MindMeld24 Analysis.exe for MindMeld24 AF
- or*
- navigating to the MindMeld24 folder and double-clicking on an .em3 scheme file.

Additionally, you can configure your PC to start a MindMeld24 application by:

- incorporating MindMeld24 on your Windows taskbar; select the MindMeld24 folder from the Programs listing, then select MindMeld24 LDC or MindMeld24 AF.
- creating shortcuts to place the MindMeld24 icons on your Windows desktop.

Once started, the selected MindMeld24 main window is displayed. If you started MindMeld24 by double-clicking a scheme file, the selected scheme is automatically loaded and displayed.

*Ensure that you have only a single copy of each MindMeld24 application running since a Windows O/S does allow multiple copies of an application to be open.*

## 2.4 SYSTEM CHECK-OUT

*These instructions are intended as a System check-out procedure only. They are not intended to convey any details of the tools or setting-up schemes.*

If you have not already done so, follow the instructions in the Mindset24 Hardware Reference Manual to prepare the hardware.

With the Mindset24 hardware connected and turned-on, and MindMeld24 running for this check-out procedure, you are accessing a pre-configured scheme consisting of two tools, the MS24 and the Stripchart. *With the scheme operating*, the Stripchart tool displays the raw data stream from the Mindset hardware much as a paper stripchart instrument does.

### Mindset/MindMeld24 Check-out Procedure Instructions

<b>User Action</b>	<b>Result/Comments</b>
Access <i>Scheme</i> → <i>Load</i> .	The Select the Scheme window displays.
Navigate to your MindMeld24 folder.	
Navigate to your Sample Scheme subdirectory.	
Select the sample Stripchart.em3 scheme.	The MindMeld24 LDC main window displays.
Access <i>Scheme</i> → <i>Start Scheme</i> or click the Start button if you have the toolbar toggled on.	<i>You should see Mindset s SCSI light illuminate. This indicates that Mindset is sampling and sending data to MindMeld24. If the SCSI light does not illuminate, check to ensure that Mindset is installed correctly and the power is on. Follow instructions in the Mindset24 Hardware Reference Manual.</i>
<i>Lightly</i> tap your fingers on Mindset s input jacks.	Observe the results of this tapping on the stripchart s traces.
Access <i>Scheme</i> → <i>Stop Scheme</i> or click the Stop button if you have the toolbar toggled on.	<i>Observe that Mindset s SCSI light goes out.</i>

Try starting and stopping the scheme several times. When you are finished, make sure the scheme is stopped.

At this point, you have started the MindMeld24 application and verified that the Mindset hardware is connected and operating correctly.

### 3 MINDMELD CONCEPTS

The following section instructs you on the basic concepts of data capture and how to use the MindMeld24 GUI.

#### 3.1 BASIC MINDMELD24 CONCEPTS

MindMeld24 is used for both live EEG data capture and recorded playback of data through a variety of visual analysis tools. These tools are selected and configured using a methodology referred to as *schemes*. A MindMeld24 scheme is the way to specify which tools you want to use, how each tool is configured and how the tools are connected to one another. Each tool is described in detail in the next section.

You can use one of the pre-configured (see Appendix C for a listing) schemes, edit them to suit your particular requirements or create new ones. Once created, a scheme can be saved as a file and reused later. You can recognize scheme files by their .em3 suffix.

MindMeld24 data are handled internally as either *streaming* or *array*.

##### 3.1.1 Streaming Data

*Streaming data are the raw data collected from the Mindset hardware. These are the type of data seen on a stripchart. Each datum is a 16 bit binary value representing the amplitude of the EEG signal at a specific point in time.*

For example, if Mindset's sampling rate is set to 256 samples per second (SPS), there are 256 of these 16 bit values collected every second from each channel. Since Mindset24 has 24 channels, the total number of samples collected in one second is  $256 \text{ SPS} \times 24 \text{ channels} = 6,144 \text{ samples}$ . (Each datum sample contains 16 bits.)

*Streaming data are fed from the MS24 tool into the:*

- *Stripchart tool* for visualization of the raw EEG signals.
- *FFT (Fast Fourier Transform) tool*. As its name implies, this tool performs an FFT on the streaming data to break the signal into its individual frequency components.

##### 3.1.2 Array Data

Once processed by the FFT tool, raw data are now *array data* since they are *handled in two dimensional arrays (frequency vs. amplitude)*. Amplitude is measured in Mindset/MindMeld from peak-to-peak in microvolts (V).

Several tools (e.g., Compressed Spectral Array, Histogram and Neuromap) require array data as their input since they perform analysis using the frequency components of the EEG signals.

*Keep these two types of data in mind when you are creating or editing your schemes.*

## 4 THE MINDMELD24 LIVE DATA CAPTURE APPLICATION

MindMeld24 Live Data Capture (LDC) is the process of obtaining the original raw data from your patient and saving it to your hard drive if required.

### 4.1 MINDMELD24 LIVE DATA CAPTURE MAIN WINDOW

These paragraphs describe how to use MindMeld24 LDC s main window from the perspective of the GUI (menu options, toolbar and window controls).

#### 4.1.1 MindMeld24 Live Data Capture Main Window Controls

The MindMeld24 LDC main window contains a main menu, a toolbar and a display area. Buttons in the *toolbar* are shortcuts for selecting the most frequently used menu options. The *display area* is where scheme tools are placed.



The MindMeld24 LDC Main Window

## **File Menu**

### **File → Exit**

Quits the MindMeld24 LDC application. You are prompted to save any open files.

## **Scheme Menu**

### **Scheme → Start Scheme**

Starts a scheme. MindMeld24 acquires live data from the Mindset hardware for presentation in the selected tools. *If you configured the Storage tool (refer to paragraph 4.3.2 for additional details) to save the raw EEG data for later replay, the specified file is opened when the scheme is started.*

### **Scheme → Pause Scheme**

Pauses current scheme operation. If the scheme is already paused, restarts the scheme. When you pause and restart a scheme, new data are appended to the end of the existing data file. *If you configured the Storage tool (refer to paragraph 4.3.2 for additional details) to save the raw EEG data for later replay, the specified file remains open when the scheme is paused.*

### **Scheme → Stop Scheme**

Stops the scheme. If you configured the Storage tool (refer to paragraph 4.3.2 for additional details) to save the raw EEG data for later replay, the specified file is closed when the scheme is stopped. *Each time you stop a scheme, the corresponding EEG data file is closed. If you subsequently restart the scheme, a new .bin file is created and the base file name is appended with a sequential number. For example; filename1.bin, filename2.bin, etc.*

### **Scheme → Edit Scheme**

Presents the Create/Edit Scheme window to create a new scheme or to modify a selected scheme. Refer to paragraph 5.1 for detailed information.

### **Scheme → Show Tools**


Displays the tools from the currently selected scheme in the display area. Use this menu option to redisplay hidden tools. (Tools are discussed in paragraphs 4.2 and 4.3.)

### **Scheme → Hide Tools**

Removes all tool windows from the display area.

### **Scheme → Load Scheme**


Opens a standard Windows file selection window titled `Select the Scheme file` to select previously created and saved schemes. Schemes have the .em3 file extension. Select the scheme to load and click Open. The selected scheme displays with its appropriate tools. Any previously loaded scheme is closed.



Select the Scheme File Window

### Scheme → Save Scheme

Opens a standard Window O/S file selection window titled *Save the Scheme file* to save the scheme currently under development. *The scheme must first be Accepted in the Create/Edit Scheme window* (refer to paragraph 5.1). Enter a name in the *File name:* field. The .em3 file extension is automatically appended when you click Save.



Save the Scheme File Window

## View Menu

### View → Toggle the Toolbar

Use this option to show or hide (i.e., *toggle*) the toolbar.

The *toolbar* contains the following buttons which correspond to several Scheme menu options described in detail above. Some buttons are grayed-out until their function is appropriate.

- Load Scheme: displays the Select the Scheme File window to select the scheme to open.
- Save Scheme: displays the Save the Scheme File window to save an Accepted scheme.
- Start Scheme: starts the scheme (i.e., data acquisition).
- Pause Scheme: pauses or restarts the scheme (i.e., data acquisition).
- Stop Scheme: stops the scheme (i.e., data acquisition).

## Special Menu

### Special → Calibrate

This option refers to the optional Calibrator hardware. When your Mindset hardware is calibrated, scalars (multipliers) are determined for each channel and saved in your PC's System Registry. Refer to Appendix B for additional details.


### Special → Use Calibration Scalars

Applies the calibration scalar values to incoming raw data prior to displaying the data in any tool window(s). Refer to Appendix B for additional details.

### Special → Audio Panel

The Audio Stimuli Control Panel provides audio stimuli to a patient while an EEG is in progress. *To use this feature you must have a sound card capable of pulse code modulation (PCM) audio.* Most sound cards support this capability.

It is best to provide audio stimulation using a Nolan Computer Systems, L.L.C. (NCS) Echofone product (details are at [www.echofone.com](http://www.echofone.com)). The piezo-electrodes on the Echofone device permit precise placement for high frequency injection, beat modulation and bone conduction. You can also use regular stereo headphones (plugged directly into your sound card) placed on your patient.



Audio Stimuli Control Panel

The left and right channel panes control the frequency and volume of the sound sine-wave. A channel is on if a checkmark appears in the "Channel On" box. To turn the channel off, click in the "Channel On" check box to remove the checkmark. To turn the channel back on, click in the check box again.

- Enter a frequency from 20 Hz to 5000 Hz.
- Adjust each channel's volume using the slide control.

To start the audio, click the "Turn Audio On" button. This button acts as a toggle to turn audio the on and off. When audio is on, the text in the button changes to "Turn Audio Off".

## Help Menu

### Help → About MindMeld24

Displays information about MindMeld24 LDC, including the software version number. The About window also has a button to *toggle on or off the background MindMeld24 bitmap shown on the main window.*

### Other MindMeld24 LDC Main Window Capabilities

In the MindMeld24 LDC main window, *right-click on any displayed tool window to show a pop-up menu.*

- Hide Window: removes that individual tool window. (All displayed tool windows are removed simultaneously using the *Scheme → Hide* option.)
- Edit Scheme: displays the Create/Edit Scheme window.

Tool windows that display streaming and array data can be resized. Other windows; such as the Storage Tool window, are fixed in size.

Once a scheme is established, each individual tool window can be moved around within the MindMeld LDC main window. Move/resize the tool window(s) into position before you start a scheme. Tool window size(s) and position(s) are preserved when a scheme is saved.

## 4.2 MINDMELD24 LIVE DATA CAPTURE TOOLS OVERVIEW

The remainder of this section describes each MindMeld24 LDC tool. These tools interact with each other and with the Mindset hardware to provide various types of visualization and analysis information. What follows first is a short description of each tool and a graphic of its icon is shown. Then each tool and its properties is explained in detail. Since the tools interact, it is necessary to have a basic familiarity with all the tools before providing an in-depth presentation of each tool.

The tools are discussed in the order they appear in the Create/Edit Scheme window (refer to paragraph 5.1).

### Mindset MS24 Tool Overview




The Mindset MS24 acquisition tool *communicates with and controls the Mindset24 hardware.* Selections within this tool set hardware parameters such as Mindset s SCSI identification (ID) number and the sampling rate. In this tool you also specify montage channel label assignments.

## Storage Tool Overview



The Storage Tool saves live EEG data to a file on your hard disk drive. Live EEG data is not saved for later replay and analysis unless a *named save file* is specified in this tool.

## Stripchart Tool Overview



The Stripchart tool *presents an EEG data stream just like an old fashioned paper stripchart machine*. You can configure the Stripchart tool to display any or all twenty-four channels, change the displayed amplitude and set the horizontal time scale. You can also specify the arrangement of the displayed channels.

## Fast Fourier Transform (FFT) Tool Overview



This tool *accepts a raw data stream as input*, performs an FFT on the data and *outputs a data array* (a signal which has been decomposed into its frequency and amplitude components.)

## Averager Tool Overview



The Averager tool *accepts array data* as input and averages these data together over a specified number of spectra. This is handy for filtering out brief transients in the data in order to visualize longer-term trends.

## Histogram Tool Overview



The Histogram tool *accepts array data* as input and outputs these data to a histogram. Use a histogram to visualize the frequency and amplitude components of the EEG signal.

## Neuromapper Tool Overview




The Neuromapper tool *accepts array data* as input and outputs these data to a neuromap, a color-coded display of brainwave activity by frequency and amplitude superimposed onto a stylized drawing of a patient's head.

## Compressed Spectral Array (CSA) Tool Overview



The CSA tool *accepts array data* as input and outputs these data to a CSA Tool window to observe how the frequency and amplitude components of an EEG signal develop over time. This is handy for visually spotting trends in the data.

## Hemispheric Tool Overview



The Hemispheric tool *accepts array data* as input and outputs these data to a hemispheric display which graphically shows brain wave activity by band on each brain hemisphere.

### 4.3 TOOLS IN DETAIL

Each tool has a primary properties window in which the tool is configured. Several tools have additional subwindows for setting other tool parameters. All of the configuration options are discussed in the paragraphs below. Following each tool's properties discussion, its corresponding window is covered.

*Right-clicking on an activated tool icon in the Create/Edit Scheme window displays that tool's properties window.* Refer to paragraph 5.1 for additional details.

MindMeld constrains the data entry fields so that you cannot enter improbable values.

#### Data Presentation


- If you invoked the *Special → Use Calibration Scalers* to the incoming raw data, scaled data is shown in the appropriate tool window(s).
- If you did not invoke the *Special → Use Calibration Scalers* to the incoming raw data, un-scaled data is shown in the appropriate tool window(s).

Additional details on calibrated data are discussed in Appendix B.

#### 4.3.1 Mindset MS24 Tool

The Mindset MS24 data acquisition tool is used *to communicate with and control the Mindset24 hardware*. Data are collected from the hardware and made ready to interact with the other tools.

*This tool is automatically activated in the Create/Edit Scheme window since the Mindset hardware is a required component for schemes to function.* Right-clicking the MS24 tool icon displays the MS24 Tool Properties window.



The MS24 Tool Properties Window

**Host Adapter** pane: set the host adapter number for your Mindset. Also review the Auto Scan button discussion below.

**SCSI ID** pane: when you installed your Mindset hardware, you specified a unique SCSI identification (ID) number. Ensure that the number in this pane matches the SCSI ID number of your Mindset hardware. If the numbers are not the same, MindMeld is unable to communicate with Mindset. Also review the Auto Scan button discussion below.

**Auto Scan** button: click the Auto Scan button to have MindMeld search for the Mindset hardware. *If the Mindset hardware is found, the host adapter and SCSI ID numbers are set with the appropriate numbers.*

**Sample Rate** pane: Mindset hardware is able to sample at 64, 128, 256, and 512 samples per second (SPS). (Refer to paragraph 3.1.1 for additional details.) A sample rate of 256 SPS is sufficient for most EEG requirements. Use higher sample rates to more accurately capture the patient's data. *Note that a high sample rate results in a very large data file.*

**Block Size** pane: the default value is 768 bytes per block. The Block Size setting is provided to help eliminate SCSI bus noise which may intrude onto the EEG data. *You should not have to reset the Block Size unless you suspect that SCSI bus noise is showing up in the data.*

With some SCSI cards the noise is more pronounced. In general, if a SCSI card draws its termination power from Mindset, the noise is more likely to show up. *It is always best to connect Mindset as the only device on your SCSI bus chain.* The SCSI card in your PC should be supplying the termination power. If this is your setup, you can turn Mindset's termination power off and remove the termination block (refer to the Mindset24 Hardware Reference Manual). These steps may eliminate the noise problem.

Depending upon several factors (such as sample rate, SCSI bus termination, SCSI bus configuration, etc.), it is possible for SCSI bus noise to intrude upon EEG data. This noise usually shows up as a spurious 2, 4, 8, 16 or 32 Hertz (Hz) signal. This is because MindMeld may be transferring data from Mindset at these intervals and the SCSI bus then has heavy activity. The SCSI bus noise is best seen by taking a minute or two of data without any other connections (i.e., an electro-cap, etc.) to Mindset's front panel. If you observe such noise, try setting the Block Size to a smaller value. This forces MindMeld to access Mindset at a higher rate because the same amount of data is transferred each second. The idea is to make this access rate higher than 32 Hz, thus moving it out of Mindset's passband. For example, assume that you are sampling at 256 SPS and the SCSI block transfer size is set to 768 bytes. Every second, there are 12,288 bytes of data transferred across the SCSI bus:


$$256 \text{ samples/second} \times 2 \text{ bytes/sample} \times 24 \text{ channels} = 12,288 \text{ bytes/second.}$$

Since MindMeld captures these data in 768 byte blocks (the SCSI Block Size), the access frequency from MindMeld to Mindset is  $12,288 / 768 = 16$  accesses per second. Under some circumstances this may show up as a spurious 16 Hz signal in the EEG data. In the above example, if you were to reset the Block Size to 384 bytes from 768, the access rate would be forced to double to 32 accesses per second, thus moving the noise signal to the upper end of the passband. Refer to the Mindset24 Hardware Reference Manual for further discussion.

**Channel Label Assignment** pane: channel label assignment mapping depicts which Mindset input channel is connected to which montage point of the Electro-cap. Some of the tools use the labels entered here as default labels on their windows. *This version of MindMeld assumes the default linked-ear montage for all neuromapping.*

All the information set in the Mindset MS24 Tool Properties window is saved when you click OK to save a scheme.

When any scheme is invoked, the Mindset MS24 Tool window automatically displays.



Mindset MS24 Tool Window

**Mindset properties** pane: your selected properties are listed.

**Mindset Internal Buffer (percent full)** bar: this progress bar fills to alert you to Mindset s internal buffer usage. The progress bar is:

- *green* when the buffer is less than 59% full
- *yellow* between 60 - 79%
- *red* at 80%.

*If the buffer overflows, the scheme stops and an error dialog box alerts you to use fewer tools in the scheme or a slower sample rate. This is not likely occur on PCs faster than 400 MHz.*

**Run time:** provides the elapsed time that this scheme has been running. Pausing the scheme continues to run the clock after the scheme is restarted; while stopping the scheme resets the clock.

### **4.3.2 Storage Tool**

The Storage tool saves live EEG data to a file. Live EEG data is not saved to your hard disk drive for later replay and analysis unless a *named save file* is specified in the Storage tool. To save live incoming data to a file, specify a file name and file/folder location. *Data are saved into the named file when you subsequently start the scheme.*

- If you invoked the *Special → Use Calibration Scalers* to the incoming raw data, the scaled data is saved using the Storage tool.
- If you did not invoke the *Special → Use Calibration Scalers* to the incoming raw data, the un-scaled data is saved using the Storage tool.

#### Data File Formats

MindMeld saves data in the *.bin* file format. A *.bin* extension is automatically appended to your file name.

## Starting, Pausing, Stopping a Scheme


Data are saved into the named file on your hard disk drive when you *start* a scheme.

If you *pause and restart* a scheme, EEG data continues to accumulate in the same .bin file.

The specified file is closed when the scheme is stopped. *Each time you stop a scheme, the EEG data file is closed. When you subsequently restart the scheme, a new .bin or file is created and the base file name is appended with a sequential number.* For example; filename1.bin, filename2.bin, etc. The initial file remains intact and accessible for replay.


(Return to paragraph 4.1.1 for information on starting, pausing, stopping and restarting schemes.)

Right-clicking the Storage tool icon in the Create/Edit Scheme window displays the standard Windows O/S save file window. For MindMeld purposes, this window is titled *Save the raw EEG data (suffix will be added automatically)* and represents the Storage Tool Properties window.



Storage Tool Properties Window

When any scheme is invoked for which you indicated on-line storage, the Storage Tool window automatically displays.



Storage Tool Window

**File name:** shows the previously entered file name for this scheme.


**Samples written:** the number of data samples that this scheme has accumulated. For the Mindset24, each data sample is 48 bytes (*24 channels x 2 bytes/channel*).

**Bytes written:** the total number of bytes accumulated in the file since the scheme has been running.

### 4.3.3 Stripchart Tool

The Stripchart tool presents a raw EEG data stream just like a paper stripchart machine. The EEG data is live, raw data during a data capture session *or* file data during an off-line analysis session. *Stripchart data are presented in real-time* during a live data capture session. For this reason, in the Create/Edit Scheme window discussed later, you can only connect a Stripchart tool directly from the MS24 tool to receive streaming data, you cannot connect it to any other tool.

The Stripchart Tool Properties window is used to set the attributes of the display when the Stripchart tool is part of a scheme. You can have two Stripchart tools running at one time, each with a different configuration. Stripchart tool(s) are invoked in the Create/Edit Scheme window by right-clicking on the Stripchart tool icon.



Stripchart Tool Properties Window


**Number of Channels:** set the number of channels to display. *The default is all 24 channels.* If you display fewer than 24 channels, all trace data still accumulates to the save file named in the Storage tool, if invoked.

**Sweep Time (sec):** set the amount of time for the traces to sweep across the entire window. For example, setting the sweep to 16 accumulates 16 seconds of data on the screen. *The default is 8 seconds.*

**Channel Properties** pane: set unique information about each channel or make all of the channels the same configuration.

- Channel Number: select the desired channel number by accessing the up/down arrows.
- Channel Label: enter your preferred label text for this channel number.
- Amplitude ( V/pix): set the displayed amplitude in microvolts per pixel ( V/pix). *This number represents the number of microvolts required to cause the traces to vary by one vertical pixel.* For example, setting 5 V/pix means that 1 vertical screen pixel represents 5 microvolts of signal.

- ◆ The *lower* the V/pix number, the larger the displayed signal amplitude; that is, the *display is more sensitive*.
- ◆ The *higher* the V/pix number, the lower the displayed amplitude; that is, the *display is less sensitive*.
- Trace Color: click the Edit button to access the standard Windows O/S Color window to select the color for the channel s tracing on the displayed stripchart.




Windows O/S Color Window

- Set All Channels: clicking this button sets all 24 channels to the same parameters (i.e., amplitude and trace color) as the selected channel.

**Stripchart Name:** assign a name to this stripchart in the data entry field.

**Channel Labels** pane: these options effect the display of the channel labels (refer to the Channel Properties pane discussion above) on the stripchart. *Be aware that certain color combinations for the label text and background color are visually more readable than others.*

- Font ANSI\_CHARSET: click the Change button to access the standard Windows O/S Font window to select the font for the channel label. The American National Standards Institute character set (ANSI\_CHARSET) is the default font.




Windows O/S Font Selection Window

- Ch. Label Color: click the Edit button to access the standard Windows O/S Color window to select the *color for the channel label s text*.

- Ch. Bkgd. Color: click the Edit button to access the standard Windows O/S Color window to select the *background color for the channel label*.

**Chart Appearance** pane: these options effect the appearance of the trace data on the stripchart. *Be aware that certain color combinations of trace, ruler and background color are visually more readable than others.* (Refer to the Channel Properties pane discussion above for setting color information.)

- Trace Bkgd. Color: click the Edit button to access the standard Windows O/S Color window to select the *color for the background on the stripchart*. (When printing a stripchart you may want to have a white background to save time and ink or toner.)
- Vertical Rule Color: click the Edit button to access the standard Windows O/S Color window to select the *color for the vertical rule* displayed on the stripchart.
- Horiz. Rule Color: click the Edit button to access the standard Windows O/S Color window to select the *color for the horizontal rule* displayed on the stripchart.
- Show Channel Labels: click in the check box to show channel labels on the displayed stripchart. Uncheck this box remove channel labels.
- Show Vertical Rule: *each vertical rule represents one second of time*. Click in the check box to show vertical rules on the stripchart. Use the pull-down option to select the rule style. Uncheck this box to remove vertical rules.




Stripchart Tool Properties Chart Appearance Pane

- Show Horizontal Rule: click in the check box to show horizontal rules on the stripchart. Use the pull-down option to select the rule style. Uncheck this box to remove horizontal rules.

**Assign Channel Map** : displays the Stripchart Channel Map window. In this window, select the channel s trace display position to correspond with each hardware channel for the stripchart. This becomes useful if you are observing fewer than 24 channels on the stripchart. *Clicking on each numbered Trace button* displays the Assign Position window in which to select any hardware channel to show at the selected trace display position on the stripchart. You may set all traces to the same hardware channel.




Stripchart Channel Map Window



Stripchart Trace Assign Position Window

All the information set in the Stripchart Tool Properties window is saved when you save a scheme.

When any scheme is invoked for which you indicated a stripchart(s), Stripchart tool window(s) automatically display with the appropriate properties. The window title indicates which slot position(s) (i.e., Stripchart 1 and Stripchart 2) are being displayed. Refer to paragraph 5.1 for a discussion on reference slots.



Stripchart Tool Window

#### 4.3.4 Fast Fourier Transform (FFT) Tool

The FFT tool is used to break an incoming signal into its individual frequencies that make up that signal. This tool *accepts a raw EEG data stream as input and outputs a data array*. This array is two dimensional, made up of specific frequencies and the amplitude of those frequencies. Several of MindMeld's tools use this *near real-time array data* as their input.

A rigorous discussion of FFT analysis and usage is beyond the scope of this manual. The following discussion is somewhat simplified and is intended to help you select from among the FFT tool parameters. For further details, see the FFT Properties 3.5 application program included on the MindMeld distribution disc.

##### FFT Epoch Length

This is the number of samples that the FFT uses in performing the transform. It is in effect the time record for the FFT. In general, a longer epoch results in a more accurate FFT output (especially at lower frequencies) but at the cost of increased processing time. *MindMeld uses a 128 sample epoch with a 0.5 Hz resolution.*

##### FFT Overlap

The overlap percentage specifies how much the epoch (or data blocks) are to be run together. The Fourier transform prefers the data be in one large continuous stream. In reality, this stream must be broken up into blocks (epochs) of finite size. This introduces adverse spectral effects; such as, phase errors, ringing and spectral leakage. To help overcome these effects, techniques referred to as windowing (refer to the next set of paragraphs) and block overlap are used. Overlapping the blocks helps the FFT average out certain effects and increases the rate at which the display can be updated.

For example, assume the time record is at two second intervals. MindMeld would have to wait for the entire time record to be captured before performing the FFT. The display would then be updated at two second intervals. Using overlap, data from the previous block are used to compute the FFT and MindMeld does not have to wait for the entire time record before transforming it. Hence, the display is updated faster. This is true up to a certain point, depending on the processing power of your PC.

A higher overlap percentage results in more transforms having to be performed each second and your PC may not be fast enough to do this while handling all of the other processing tasks. *MindMeld uses an overlap percentage of 75%.*

## FFT Windows

*Think of an FFT window as a filter* through which data are passed prior to computing the FFT. FFT windows are used to reduce certain spectral effects caused by performing the FFT on finite data sets. *The incoming signal is likely not be periodic within the time record, thus causing phase errors.* Also, if a signal's frequency is not precisely on a frequency bin (0.5 Hz for Mindset/MindMeld), the signal may smear across two or more frequency bins. An FFT window tapers the data to zero at each endpoint which forces it to be periodic within the time record. The different FFT window shapes determine just how the edges are tapered (steeply or gradually).

In the following discussion, the term "*frequency bin*" refers to each of the 128 frequencies (0.5 Hz to 64 Hz at 0.5 Hz intervals) computed by the MindMeld FFT. The term *selectivity* refers to the ability of the FFT to distinguish signals which are close together in frequency.

- Rectangular

This is not an FFT window (i.e., filter) at all. *The FFT is computed on the raw data.* With the rectangular window, a signal that is exactly on a frequency bin shows as a single line in the histogram tool because it is exactly periodic within the time record. *Only those frequencies which are exactly on the frequency bins show with the proper amplitude and frequency in the histogram (and other tools which accept their input from the FFT tool).* If a frequency is between bins, it smears out and affects every bin of the spectrum. For this reason, the Rectangular window is not a good option for general EEG purposes.

Since the Mindset Calibrator injects a precise 16 Hz signal into the Mindset, the special calibration scheme contains an FFT tool using the Rectangular window option.

- Hanning

The Hanning window has an amplitude variation of about 1.5 dB (decibel) for signals between bins and provides fairly good selectivity. *The Hanning window limits the performance of the FFT when two signals are very close together in frequency but different in amplitude.*


- Blackman

*The Blackman window is a very good window to use with the Mindset.* It has better accuracy than the Hanning, very good selectivity and fast filter rolloff. Signals close together can be resolved even when their amplitudes differ to a large degree.

- Flattop

The Flattop window improves on the amplitude accuracy of the Hanning window. The Flattop window's between bin amplitude variation is about .02 dB but its selectivity is a bit less precise. The Flattop window has a wide passband and very steep rolloff characteristics. *Therefore, signals appear wide but do not leak across the entire spectrum.*

Use the FFT Tool Properties window to set the FFT tool's parameters. FFT tool(s) are invoked in the Create/Edit Scheme window by right-clicking on the FFT tool icon.




FFT Tool Properties Window

**FFT Window:** select the type: Rectangular, Hanning, Blackman or Flattop from the pull-down menu option.

**FFT Window scale factors** pane: enter the scale (i.e., weighting) values you want. *These are multipliers used to bring the resulting amplitude signals into compliance with other databases that you may be using.* Also, since the different FFT window types result in a reduction of the actual signal amplitude, these factors are used to bring the displayed amplitude up to the true amplitude. *The values shown in the above figure are calibrated values determined by injecting a signal of known frequency and amplitude into the Mindset and observing the FFT output.*

All the information set in the FFT Tool Properties window is saved when you save a scheme.

When any scheme is invoked for which you indicated an FFT, the FFT Tool window automatically displays with your selected properties listed. The window title indicates which slot position(s) (i.e., FFT 1 and FFT 2) are being displayed. Refer to paragraph 5.1 for a discussion on reference slots.



FFT Tool Window


### 4.3.5 Averager Tool

The Averager tool *takes FFT array data as input and averages these data over time. Brief transients are smoothed out.* This facilitates visualizing long-term trends in the data.

Because of this averaging, *the data display within tool window(s) connected to the Averager tool are initially delayed and slower in comparison to a Stripchart tool which takes data directly in real-time.*

Outputs from the FFT to be averaged together before being displayed are selectable from 2 to 16 spectra. *A spectra is one output of a sample block (i.e., 128 samples in length) from the FFT.* The higher the number of spectra to be averaged, the smoother any rapid transitions (e.g., artifacts) display on-screen.


Averager parameters are set in the Averager Tool Properties window. Averager tool(s) are invoked in the Create/Edit Scheme window by right-clicking on the Averager tool icon.



Averager Tool Properties Window

All the information set in the Averager Tool Properties window is saved when you save a scheme.

When any scheme is invoked for which you indicated an Averager, the Averager Tool window automatically displays with your spectra selection shown. The window title indicates which slot position(s) (i.e., Averager 1 and Averager 2) are being displayed. Refer to paragraph 5.1 for a discussion on reference slots.




Averager Tool Window

### 4.3.6 Histogram Tool

The Histogram tool permits visualization of spectral EEG data (the frequency and amplitude components of the EEG signal). It *accepts FFT array data from the FFT tool or Averages tool and outputs these data to a histogram.*

Histogram parameters are set in the Histogram Tool Properties window. Histogram tool(s) are invoked in the Create/Edit Scheme window by right-clicking on the Histogram tool icon.



Histogram Tool Properties Window

**Number of channels** pane: set the number of channels to be displayed. *If you display less than 24 channels, all data still accumulates to the save file named in the Storage tool, if invoked.*

**Display band** pane: determines the frequencies that are displayed. You can limit the amount of data displayed by selecting a narrower range of frequencies. *Each bar in a histogram represents the amplitude of signals which fall within that band.* For example, if you select the Passband option, the amplitude of the bar varies in accordance with frequencies in the 2 Hz to 34 Hz range.


Select fixed parameters (Passband or Full scale) for the range of frequencies to be displayed or customize the range by selecting Other.

- Passband (2 Hz to 34 Hz): click to display the listed frequencies.
- Full scale (1 Hz to 63 Hz): click to display the listed frequencies.
- Other: click to specify your custom range between 1 Hz and 63 Hz:
  - ◆ Low freq: enter the lowest frequency to be displayed.
  - ◆ High freq: enter the highest frequency to be displayed.

**X-Axis Labels** pane: select whether the X-axis labels are displayed.

**Maximum amplitude** pane: enter the maximum vertical amplitude that can be displayed. The amplitude may be from 1 V to 120 V.


**Channel Map** pane: select the channel s display position to correspond with each hardware channel for the displayed histogram. *Clicking on each numbered Position button displays the Assign Position window* in which to select the hardware channel to show at the selected numbered display position.



Assign Position Window for Histogram Tool

All the information set in the Histogram Tool Properties Window is saved when you save a scheme.

When any scheme is invoked for which you indicated a histogram, the histogram automatically displays with the appropriate properties shown. The Histogram Tool window title identifies the reference slot (refer to paragraph 5.1) and maximum voltage (from your Maximum amplitude selection above) for this histogram. *The vertical scale is in microvolts with 10 divisions marked.* Histogram colors are fixed for each band.



Histogram Tool Window (4 Channel Example)

### 4.3.7 Neuromapper Tool

The Neuromapper tool enhances *visualization of spectral EEG data as a neuromap*. Neuromaps are helpful for visualizing localized EEG information. The various EEG waveforms are color-coded and shown at their respective electrode sites on a stylized drawing of a patient's head. *You can have several Neuromap tool windows open simultaneously, each displaying a different frequency band.*

The neuromapper tool is currently hardwired to specific channels and electrode positions and assumes the standard 10-20 configuration. Refer to Chapter 2, the electrode configuration paragraphs of the Mindset24 Hardware Reference Manual for details.

The Neuromapper tool *accepts array data from the FFT tool or Averages tool* and outputs these data to a composite neuromap. That is, all frequency bands and amplitudes are shown on the same neuromap.

Neuromap parameters are set in the Neuromap Tool Properties window. Neuromap tool(s) are invoked in the Create/Edit Scheme window by right-clicking on the Neuromap tool icon.



Neuromap Tool Properties Window

**Threshold** pane: select the threshold for the brightness of the display. *This amplitude represents the lowest signal amplitude which displays. Any input amplitudes below this level display in black.*

**Range** pane: set the range (sensitivity) of the neuromap. *A setting of 20 V means that it takes 20 V of signal (above the threshold setting) to illuminate the neuromap at the brightest colors.* The Range setting is added to the Threshold setting to determine the overall amplitude range that is displayed.


For example, selecting a Threshold of 10 V and a Range of 80 V means that signal amplitudes less than 10 V are shown as black and signal amplitudes of 90 V and above are shown at full illumination.

All the information set in the Neuromap Tool Properties window is saved when you save a scheme.

When any scheme is invoked for which you indicated a neuromap, a neuromap automatically displays with the appropriate properties. The Neuromap Tool window title identifies its reference slot (refer to paragraph 5.1). The legends are:

- Left-side: display scale in microvolts based on your tool property settings.
- Right-side: band colors are pre-set to the following:

<i>Band Frequency (Hz)</i>	<i>Color</i>
Delta 1.0 — 3.0	Blue
Theta 3.5 - 7.5	Green
Alpha 1 8.0 — 14.0	Yellow
Alpha 2 14.5 — 18.0	Orange
Beta 18.5 — 40.0	Red




Neuromap Tool Window

### 4.3.8 Compressed Spectral Array (CSA) Tool

The Compressed Spectral Array (CSA) tool permits *visualization of spectral EEG data over time*. The data are displayed with *frequency* along the horizontal (X) axis, *amplitude* in microvolts along the vertical (Y) axis and *time* along the diagonal (Z) axis. (Refer to the CSA Tool window figure below.)

The CSA tool *accepts array data from the FFT tool or Averager tool* and outputs these data to a CSA tool window.

CSA parameters are set in the CSA Tool Properties window. CSA tool(s) are invoked in the Create/Edit Scheme window by right-clicking on the CSA tool icon.



CSA Tool Properties Window

**Number of channels** pane: set the number of channels to be displayed. *If you display less than 24 channels, all data still accumulates to the save file named in the Storage tool, if invoked.*

**Display band** pane: determines the frequencies that are displayed. You can limit the amount of data displayed by selecting a narrower range of frequencies. Each trace in a CSA tool window represents the amplitude of signals which fall within that band. For example, if you select the Passband option, the amplitude of the trace varies in accordance with frequencies in the 2 Hz to 34 Hz range.


Select fixed parameters (Passband or Full scale) for the range of frequencies to be displayed or customize the range by selecting Other.

- Passband (2 Hz to 34 Hz): click to display the listed frequencies.
- Full scale (1 Hz to 63 Hz): click to display the listed frequencies.
- Other: click to specify your custom range between 1 Hz and 63 Hz:
  - ◆ Low freq: enter the lowest frequency to be displayed.
  - ◆ High freq: enter the highest frequency to be displayed.

**Maximum amplitude** pane: enter the maximum vertical amplitude that can be displayed. The amplitude may be from 1 V to 120 V.

**Colors** pane: click the Edit button to access the standard Windows O/S Color window (refer to the Channel Properties pane discussion in paragraph 4.3.3) to select the background, foreground and pen colors for a CSA tool. *Be aware that certain color combinations for the foreground and background color are visually more readable than others.*


**Channel map** pane: select the display position to correspond with each hardware channel. *Clicking on each numbered Position button displays the Assign Position window in which to select the hardware channel to show at the selected numbered display position.*



Assign Position Window for CSA Tool

All the information set in the CSA Tool Properties window is saved when you save a scheme.


When any scheme is invoked for which you indicated a CSA, the CSA tool window automatically displays with the appropriate properties shown. The CSA tool window title identifies the reference slot (refer to paragraph 5.1), maximum amplitude (i.e., the vertical scale) for the scale and frequency band (i.e., the horizontal scale) information for this CSA.



CSA Tool Window (4 Channel Example)

If you select a Histogram tool and a CSA tool, and place these windows side-by-side on-screen, you can see the frequency and amplitude of the histogram reflected as the front-most time-line within the CSA Tool window.

The CSA tool presents the best visualization of the difference between real-time FFT tool array data (refer to paragraph 3.1.1) and Averager tool near-real-time data (refer to paragraph 3.1.2) input. Although the graphic below is small, you should be able to *detect the slower presentation and smoothing effects of the Averager tool on the CSA tool display.*




Comparison of FFT Tool and Averager Tool Effects on a CSA Tool

### 4.3.9 Hemispheric Tool

The Hemispheric tool presents a stylized version of brain wave activity by band on each brain hemisphere. The Hemispheric tool *accepts array data* as input and outputs these data to a hemispheric display. There are 16 channels displayed (8 channels on each side).

Hemispheric parameters are set in the Hemispheric Tool Properties window. Hemispheric tool(s) are invoked in the Create/Edit Scheme window by right-clicking on the Hemispheric tool icon.



Hemispheric Tool Properties Window

**Maximum amplitude** pane: enter the maximum vertical amplitude from 1 V to 120 V to displayed.

**Display band** pane: determines the range of frequencies that are displayed. You can limit the amount of data displayed by selecting a narrower range of frequencies. Each bar in a hemispheric display represents the amplitude of signals which fall within that band. For example, if you set 1 for the low frequency and 2 for the high frequency, the amplitude of the bar varies in accordance with frequencies in the 1 Hz to 2 Hz range.

Select fixed parameters (Passband or Full scale) for the range of frequencies to be displayed or customize the range by selecting Other.

- Passband (2 Hz to 34 Hz): click to display the listed frequencies.
- Full scale (1 Hz to 63 Hz): click to display the listed frequencies.
- Other: click to specify your custom range between 1 Hz and 63 Hz:
  - Low freq: enter the lowest frequency to be displayed.
  - High freq: enter the highest frequency to be displayed.


All the information set in the Hemispheric Tool Properties window is saved when you save a scheme.

**Channel map** pane: determines which channel displays in which position. Default is hemispheric.

When any scheme is invoked for which you indicated a hemispheric presentation, a Hemispheric tool window automatically displays with the appropriate properties shown. The Hemispheric Tool window title identifies the reference slot (refer to paragraph 5.1) and frequency band information for this presentation. (The color legend is the same as for a Neuromapper tool window.) Band colors are pre-set to the following:

<i>Band Frequency (Hz)</i>	<i>Color</i>
Delta 1.0 — 3.0	Blue
Theta 3.5 - 7.5	Green
Alpha 1 8.0 — 14.0	Yellow
Alpha 2 14.5 — 18.0	Orange
Beta 18.5 — 40.0	Red

The electrode labels/positions (set in the MS24 Tool Properties window; refer to paragraph 4.3.1) are shown along the vertical sides with the maximum amplitude shown on the bottom axis. The vertical centerline represents the division between the right- and left-brain hemispheres.



Hemispheric Tool Window

## 5 CREATE/EDIT/RUN A MINDMELD LIVE DATA CAPTURE SCHEME


Creating and editing a MindMeld24 LDC scheme consists of linking the various tools together to display your acquired EEG data. This section explains the method used to select and connect tools into a scheme. A scheme may have from two to all of the tools within it.

Running a scheme is the process of live data capture of a patient's EEG information.

### 5.1 USING THE CREATE/EDIT SCHEME WINDOW

*Schemes are temporary until they are saved.* You can run an unsaved scheme, but if you make any changes, the scheme must be re-established the next time you want to use it.

Use *Scheme Menu* → *Edit Scheme* to display the Create/Edit Scheme window.



Create/Edit Scheme Window

**Click on a tool to place it into the scheme** pane: these are the icons for the 9 tools. You must always *click on a tool's icon first to select it*, then the icon shows in its appropriate location(s) in the Scheme tools: pane of the window. Tools do not have to be selected in icon order.

The *MS24 tool* is automatically selected since the Mindset24 hardware must always be controlled by the MindMeld software in any scheme.

You can select a number of different tools or select the same tool several times setting each invocation of that tool's parameters differently.

**Scheme tools** pane: some tools, like the Stripchart have assigned *slots (reference points for the scheme connections)*. Other tools, like the Neuro-mapper, fill from slots 10 — 15 as they are selected.

Tools must be connected together to indicate how data is to flow between the various tools. The *colored data flow arrows indicate which tools can be connected*. To connect the tools, left-click on the source tool (it highlights) and drag the cursor to the destination tool. When the destination tool highlights, release the mouse button.

- *Red* arrow: hardware, data storage and streaming data tool connections.
- *Blue* arrow: array data connections.
- *Green* arrow: array data from the FFT to Averager tool connection.

The *Create/Edit Scheme window has embedded controls which prevent you from making improper connections*. Once the tools are properly connected, an appropriately colored dotted-line displays the connection. If an improper connection is attempted, the color-coded dotted-line does not show.

If you try to connect a second input to any tool, an error dialog box displays.

### **Buttons:**

- **Accept:** click when your scheme is appropriately configured. *You must Accept a scheme first before you can Save it.* (Refer to paragraph 4.1.1.) If you Save a scheme before it has been Accept(ed), you are saving the scheme currently displayed in the display scheme area of the MindMeld LDC main window instead of the scheme you are creating/editing in the Create/Edit Scheme window.
- **Cancel:** *click to close the Create/Edit Scheme window.* A confirmation dialog box displays:
  - ◆ click Yes to close the Create/Edit Scheme window, nothing is saved
  - or*
  - ◆ click No to return to the Create/Edit Scheme window.
- **Start Over:** *click to restore the Create/Edit Scheme window to its unconnected state.* A confirmation dialog box displays:
  - ◆ click Yes to unconnect the scheme
  - or*
  - ◆ click No to return to the current scheme.

If you click Start Over and then click Accept, you clear the MindMeld LDC main window.

### Other MindMeld LDC Window Capabilities


Right-click on a tool icon to access its tool properties window.

If you need to delete a selected tool, Shift & right-click on a tool icon to remove it from the current scheme.

## 5.2 RUNNING A SCHEME

Once you have a scheme created or edited, the Mindset is on-line and your patient is properly prepared, click the toolbar Start Scheme icon or invoke *Scheme* → *Start Scheme* to begin the EEG data capture process. The tool windows display with live data and data is recorded to your hard drive if the Storage tool is part of that scheme's configuration.

Allow the System to acquire data as long as necessary to obtain an acceptable amount of data for analysis and diagnosis.



MindMeld LDC Main Window with Scheme Tools Displayed

## 6 THE MINDMELD24 ANALYSIS FUNCTIONS APPLICATION


The *MindMeld24 Analysis Functions (AF)* is used for whole session, off-line analysis. This process is often called *replay*. After EEG data has been captured and saved to your hard drive with MindMeld24 LDC, the tools in MindMeld24 AF are applied to the data to provide an overall interpretation of the patient's session.

When a data file is captured, it likely contains various artifacts; such as, eye and muscle movement, responses to unintended stimuli, etc. In MindMeld24 AF you can remove these artifacts from the statistical and analysis tool's calculations and displays. These artifact removal sections are called *Ignore Regions (IR)*. *The data is never removed from the actual EEG data file, it is only deleted from calculations and presentations.*

### 6.1 MINDMELD24 ANALYSIS FUNCTIONS MAIN WINDOW CONTROLS

The MindMeld AF main window contains a main menu, a toolbar and a display area.

Many of these menu and toolbar options are the same as from the MindMeld24 LDC main window. They are all repeated here, however, for continuity in this discussion.



MindMeld24 AF Main Window with Scheme Running

## File Menu


### File → Open EEG File

Presents the standard Windows O/S file selection window titled *Select the EEG data file to open* to select a saved EEG data file for replay. *These are the data files created during a MindMeld live data capture process when the Storage tool is selected and configured as part of a scheme.* MindMeld data files have a *.bin* file extension.

(Refer to the Storage tool discussion in paragraph 4.3.2 for details on the impact of starting and stopping a scheme as it relates to file names.)

Select the *.bin* file to load and click Open.

- The *file data shows within a selected (previously saved) scheme or within a newly created (unsaved) scheme.* Any previously replayed data file is closed.  
*or*
- If you have not selected or created a scheme at this point, the EEG data is not shown until a scheme is loaded.




Select the EEG Data File To Open Window

### File → Close EEG File

If Ignore Regions (IRs) have been selected or changed during this replay and you select this option, a confirmation dialog displays asking whether to *Save the changed IR?* Click Yes to save the IR information for the next replay, click No to disregard saving the IR information.

## File → Page Setup

This option presents the standard Windows O/S Page Setup window to set your printed page parameters.



Windows O/S Page Setup Window

## File → Exit

Quits the MindMeld AF application. You are prompted to save any open files that reflect changed schemes or IRs.

## Scheme Menu

### Scheme → Start Scheme

Starts a replay scheme with the selected data file. A replay starts from the time shown in the Stripchart 1 Tool window. *Use the scrollbar in the Stripchart 1 Tool window to indicate your preferred start time.* Refer to paragraph 6.2.4 for additional details.

After your file is loaded and replayed, the Statistics window remains on-screen and does not re-set (i.e., zeroize) until you start the same or another scheme with the current or a different EEG file. The Mindset Tool 1 window remains on-screen with the selected file name shown in its File properties: pane.

### Scheme → Pause Scheme

Pauses the current replay. If the scheme is already paused, restarts the scheme.

### Scheme → Stop Scheme

Stops the replay.

### Scheme → Edit Scheme

Presents the Create/Edit Scheme window to allow a new scheme to be created or to modify the current scheme. Refer to paragraph 5.1 for detailed information.

### Scheme → Show Tools

Displays the Statistics window and other MindMeld24 AF tool(s) windows for the current scheme.


### Scheme → Hide Tools

Removes the Statistics window and the MindMeld24 AF tool(s) windows from the main window.

### Scheme → Load Scheme

Opens the standard Windows O/S file selection window titled **Select the Scheme file** to select previously saved schemes. Select the scheme to load and click **Open**. The scheme displays with the appropriate tools. Any scheme already loaded is closed.


- If you already selected a data file to replay, data shows in the tool windows.
- If you have not selected a data file to replay, access **File → Open File** now. Refer to the **File → Open File** discussion above for additional details.



Select the Scheme File Window

### Scheme → Save Scheme

Opens the standard Windows O/S file save window titled **Save the Scheme file** to save the scheme under development. *The scheme must first be Accepted in the Create/Edit Scheme window* (refer to paragraph 5.1). Enter a name in the **File name:** field, the **.em3** file extension is automatically appended when you click the **Save** button.



Save the Scheme File Window

## View Menu

Use this option to show or hide (i.e., *toggle*) the toolbar display.



MindMeld24 AF Toolbar

The *toolbar* contains the following buttons which correspond to several Scheme and Special menu options described in paragraphs above or below. Some buttons are grayed-out until their function is appropriate.

- Open Scheme: displays the Select the Scheme File window to select the scheme to open.
- Save Scheme: displays the Save the Scheme file window to save the scheme under development.
- Start Scheme: starts the scheme (i.e., the replay).
- Pause Scheme: pauses or restarts the scheme (i.e., the replay).
- Stop Scheme: stops the scheme (i.e., the replay).
- Add to Ignore Regions: refer to Special Menu discussion below.
- Remove from Ignore Regions: refer to Special Menu discussion below.
- Respect Ignore Regions: refer to Special Menu discussion below.
- Clear Ignore Regions: refer to Special Menu discussion below.

## Special Menu

### Special → Use Calibration Scalers

This option refers to the optional Calibrator hardware. When the Mindset hardware is calibrated, scalars (multipliers) are determined for each channel and saved in your PC s System Registry. When checked, this menu option indicates that the scalars are to be used to adjust the data before it is passed to any tools.

If you already invoked *Special → Use Calibration Scalers* when you acquired the EEG data using MindMeld24 LDC, invoking this same option (*Special → Use Calibration Scalers*) in MindMeld24 AF has no effect. In other words, the data is not re-calibrated a second time. Additional details are provided in Appendix B.

### Special → Run at Full Speed

Runs the data replay at actual speed. If you are using a slower PC (200 MHz or less) and the EEG data being played back was sampled at 512 samples/second, your PC may be sluggish and not respond to input commands. If this is the case, play the file back at half-speed or slower.

### Special → Run at Half Speed

Runs the data replay at 1/2 the actual speed.

### Special → Run at Quarter Speed

Runs the data replay at 1/4 the actual speed.

### Special → Run at Eighth Speed

Runs the data replay at 1/8th the actual speed.

### **Special → Add to Ignore Regions**

Ignore Regions are selected by *left-clicking within the Stripchart 1 tool and dragging your cursor to the right or left over the data to be disregarded* in calculating Statistic window values and in tool displays. (The EEG data within an IR are not eliminated from the source .bin or .eeg data file.) As you drag your mouse over the data to be ignored, it highlights as visual feedback. (If you make a mistake and do not want to include your highlighted data in the Ignore Regions, simply click in the Stripchart 1 Tool window and the highlighted region returns to its configured background color.)

To set your selected data as an IR, invoke *Special → Add to Ignore Regions* or click the Add to Ignore Regions toolbar button. The highlighted portion shows with a red background as a visual indication that it is an IR.

Adding to an existing IR means to left-click and drag the mouse to/from a red IR to expand it, then click the Add to IR toolbar button or invoke *Special → Add to Ignore Regions* to include the new data.

As an EEG data file is replayed, *data within Ignore Regions are not included in any calculations and corresponding tool windows within that scheme.*

Ignore Region(s) are saved with the source EEG data files when you select *File → Close EEG File.*

### **Special → Remove from Ignore Regions**

Click the Remove from IR toolbar button or invoke *Special → Remove from Ignore Regions.* A small scissors icon displays, mouse to the required IR, then left-click within an existing IR. The data display returns to its original color from a red background.

### **Special → Respect Ignore Regions**

This option toggles the ability to *exclude* (i.e., *respect*) IR(s) when calculating Statistic window values and tool window displays. When the Respect Ignore Regions menu option or button on the toolbar is not showing a checkmark, data within IR(s) are included in scheme calculations.

### **Special → Clear Ignore Regions**

All IR(s) are erased from the display. The actual EEG data remains in the file.

## Help Menu


### Help → About MindMeld24

Displays information about MindMeld24 AF, including the software version number. The About window also has a button to toggle on or off the background MindMeld24 bitmap shown on the main window.

### Other MindMeld24 AF Window Capabilities

In any displayed tool window, *right-click to access a pop-up menu.*

- Hide Window: removes that tool window.
- Edit Scheme: displays the Create/Edit Scheme window.
- Print This Window: displays the standard Windows O/S Print window to print this tool window.



Windows O/S Print Window

MindMeld24 AF main and tool windows can be moved and resized in the same manner as MindMeld24 LDC windows.


## 6.2 MINDMELD24 ANALYSIS FUNCTIONS: CREATE/EDIT/REPLAY SCHEME

MindMeld24 AF capabilities for creating and editing a scheme are similar to the MindMeld24 LDC process for creating and editing a scheme. The run/replay scheme capabilities are identical between the two applications except for the use of Stripchart 1 as discussed below.

*In the following discussion, only differences between the two applications are addressed.*

### 6.2.1 MindMeld24 Analysis Functions Mindset 1 Tool Properties

Within MindMeld24 AF, the MS24 Tool Properties window only allows you to set channel label assignments. The remainder of the properties are retained from the original scheme.




MindMeld24 AF Mindset 1 Tool Properties Window

**Channel Assignment** pane: channel label assignment mapping depicts which Mindset input channel is connected to which montage point of the Electro-cap. *This version of MindMeld24 assumes the default linked-ear montage for all neuromapping.*

All the information set in the MS24 Tool Properties window is saved when you click OK to save a scheme.

When any scheme is invoked, the MS24 Tool window automatically displays. The window title also identifies the Mindset hardware model number which was used in the data capture process and the number of channels available on that hardware.



MindMeld AF Mindset Tool 1 Window

**File properties** pane: reflects information from the original scheme configuration and the selected file. (Refer to paragraph 4.3.1.)

- Disk File: None (if a data file has not yet been selected) or the name of the selected .bin file.
- File Size: the number of samples in the selected .bin file.
- Sample rate: the sample rate at which the selected .bin file was obtained.
- Channels: the number of channels in use when the selected .bin file was obtained.

**File completion** bar: the progress bar uses green fill to mark the pace of the replay.

**Run time:** shows the *elapsed time of the replay in seconds* up to 9 decimal places. The precision of this time information is 1 divided by the sample rate (SR); that is,  $1 \div \text{SR}$ . For example, if the SR is 256, then each sample takes  $1 \div 256^{\text{th}}$  of a second (or .00390625) to occur. When sample number 204 occurs, the elapsed Run time is  $204 \times .00390625$  or .796875 seconds.

**Sample:** shows the *total number of samples* displayed since the start of the replay.

## 6.2.2 MindMeld24 Analysis Functions Statistics Window

The Statistics window automatically displays in the MindMeld24 AF main window. After the first file is loaded and replayed, the Statistics window does not re-set (i.e., zeroize) itself until you start a scheme with the currently selected file. (Refer to the *File* → *Open File* discussion in paragraph 6.1.) (If you are running MindMeld on a slower CPU, there may be a slight delay before values show in the Statistics window.)

**NOTE:** In order to calculate Statistics values, *the MS-1000 tool must be connected to the FFT 1 tool in the corresponding Create/Edit Scheme window.*

The Statistics Window. Values are micro-volts peak-to-peak presents a table that contains the amplitude value from peak-to-peak at that time in the data file replay. At the bottom of the table, summary statistics for *Minimum*, *Maximum*, *Delta* and *Average* amplitude are calculated. The Delta row shows the difference between the Maximum and Minimum values.

You can start statistical calculations from any point within the EEG data file. Use the horizontal scroll bar on the Stripchart 1 tool window to scroll to the desired time position and then start the scheme. *Thus, you can calculate the statistics for any portion or all of a patient s session.* Refer to paragraph 6.2.4 below for additional details.

Channel	FP1	FP2	F7	F3	Fz	F4	F8	T3	C3	Cz	D4	T4	T5	P3	Pz	P4	T6	O1	O2	Ch20	Ch21	Ch22	Ch23	Ch24
2.0	14.0	5.9	4.0	4.0	4.0	4.0	4.0	4.0	4.0	4.0	4.0	4.0	4.0	4.0	4.0	4.0	5.1	4.0	4.0	4.0	4.0	4.0	4.0	14.0
2.5	11.5	4.6	3.1	3.1	3.1	3.1	3.1	3.1	3.1	3.1	3.1	3.1	3.1	3.1	3.1	3.1	3.9	3.1	3.1	3.1	3.1	3.1	3.1	11.5
3.0	9.6	3.8	2.6	2.6	2.6	2.6	2.6	2.6	2.6	2.6	2.6	2.6	2.6	2.6	2.6	2.6	3.3	2.6	2.6	2.6	2.6	2.6	2.6	9.6
3.5	8.1	3.2	2.2	2.2	2.2	2.2	2.2	2.2	2.2	2.2	2.2	2.2	2.2	2.2	2.2	2.2	2.8	2.2	2.2	2.2	2.2	2.2	2.2	8.1
4.0	7.2	2.8	1.9	1.9	1.9	1.9	1.9	1.9	1.9	1.9	1.9	1.9	1.9	1.9	1.9	1.9	2.4	1.9	1.9	1.9	1.9	1.9	1.9	7.2
4.5	6.3	2.5	1.7	1.7	1.7	1.7	1.7	1.7	1.7	1.7	1.7	1.7	1.7	1.7	1.7	1.7	2.1	1.7	1.7	1.7	1.7	1.7	1.7	6.3
5.0	5.7	2.2	1.5	1.5	1.5	1.5	1.5	1.5	1.5	1.5	1.5	1.5	1.5	1.5	1.5	1.5	1.9	1.5	1.5	1.5	1.5	1.5	1.5	5.7

MindMeld24 AF Statistics Window (Upper Portion)

325	0.7	0.3	0.2	0.2	0.2	0.2	0.2	0.2	0.2	0.2	0.2	0.2	0.2	0.2	0.2	0.2	0.3	0.2	0.2	0.2	0.2	0.2	0.3	0.2	0.6
330	0.7	0.3	0.2	0.2	0.2	0.2	0.2	0.2	0.2	0.2	0.2	0.2	0.2	0.2	0.2	0.2	0.3	0.2	0.2	0.2	0.2	0.2	0.2	0.2	0.6
335	0.6	0.3	0.2	0.2	0.2	0.2	0.2	0.2	0.2	0.2	0.2	0.2	0.2	0.2	0.2	0.2	0.3	0.2	0.2	0.2	0.2	0.2	0.2	0.2	0.6
340	0.6	0.3	0.2	0.2	0.2	0.2	0.2	0.2	0.2	0.2	0.2	0.2	0.2	0.2	0.2	0.2	0.3	0.2	0.2	0.2	0.2	0.2	0.2	0.2	0.6
Minimum	-18.8	-47.9	-54.0	-54.1	-54.0	-54.1	-54.1	-54.1	-54.0	-54.0	-54.1	-54.1	-54.1	-54.1	-54.1	-54.1	-50.3	-54.0	-54.1	-54.1	-54.1	-54.1	-54.0	-18.8	
Maximum	26.6	40.1	53.6	53.7	53.8	53.7	53.8	53.6	53.8	53.7	53.8	53.7	53.8	53.8	53.8	53.7	46.9	53.8	53.7	53.5	53.5	53.2	26.5		
Delta	45.5	7.8	0.4	0.4	0.3	0.3	0.4	0.3	0.4	0.3	0.4	0.3	0.2	0.3	0.3	4.4	0.2	0.4	0.5	0.6	0.6	0.8	45.3		
Average	3.9	43.9	53.8	53.8	53.8	53.8	53.8	53.8	53.8	53.8	53.8	53.8	53.8	53.8	53.8	53.8	48.0	53.7	53.8	53.8	53.8	53.8	53.8	3.9	


MindMeld24 AF Statistics Window (Lower Portion)

**NOTE:** If you are running MindMeld on a PC with a slower CPU (less than 400 MHz) and the sample rate on the replay data file is at 512, there may be a processing *delay* in the window displays and your ability to mouse and move/resize tool windows may be impaired while Statistic window values are calculated. In this instance, use *Special* → *Run at Half Speed* to minimize these processing delays.

### Other Statistics Window Capabilities

In the Statistics window, *right-click to access a pop-up menu.*

- Hide Window: removes the Statistics window.
- Edit Scheme: displays the Create/Edit Scheme window.
- Print This Window: displays the standard Windows O/S Page Setup window the first time this pop-up menu is invoked during a MindMeld24 AF session to set your printed page parameters. (That is, whenever you exit and restart MindMeld24 AF, the Windows O/S Page Setup window displays the first time you select this Print This Window pop-up menu option.) Click OK to display the standard Windows O/S Print window. Thereafter, only the Print window displays to invoke printing the Statistics window.
- Dump to File: displays a Save the statistics data window. This is a standard Windows O/S save file window to save an EEG data file, allowing you to save this statistical data as a tab de-limited text file to be opened and printed in Microsoft Word, WordPad or Notepad, etc. Statistics window data can also be saved and imported into spreadsheet or database applications for follow-on data manipulation. (EEG file data is only saved as a text file.)



MindMeld24 AF Save the Statistics Data Window

### **6.2.3 MindMeld24 Analysis Functions Storage Tool**

The Storage tool is not available when using MindMeld24 AF.

### **6.2.4 MindMeld24 Analysis Functions Stripchart Tool**

In MindMeld24 AF, the Stripchart 1 tool has special importance. *The horizontal scroll bar on the Stripchart 1 Tool window controls the position within an EEG data file from which data are played-back and Ignore Regions are set.*

Once an EEG file is open for playback, use the scroll bar to move (in one second increments) to the desired position within the file. *When the scheme is started (Scheme → Start Scheme), playback begins at the time position shown on this stripchart.*

### **6.2.5 MindMeld24 Analysis Functions Remaining Tools**

The other tools and their windows remain the same as those in the MindMeld24 LDC function, as discussed in paragraph 4.3.

# Appendix A

## Acronyms and Glossary

.bin	data file extension for new MindMeld files
.eeg	data file extension for older Mindlab/Wavelab files
.em2	MindMeld16 scheme file extension (for Mindset model MS-1000)
.em3	MindMeld24 scheme file extension (for Mindset model Mindset24)
.reg	Windows O/S RegEdit file extension
ADC	Analog-to-Digital Converter
ANSI	American National Standards Institute
CD-ROM	Compact Disc — Read-only Memory
CPU	Central Processing Unit
CSA	Compressed Spectral Array
dB	decibel
EEG	electroencephalography
FDA	Food and Drug Administration
FFT	Fast Fourier Transform
GUI	Graphical User Interface
Hz	Hertz
ID	Identification
IR	Ignore Region
MindMeld24 AF	MindMeld24 Analysis Functions application
MindMeld24 LDC	MindMeld24 Live Data Capture application
MHz	Mega-Hertz (millions of Hertz)
MS-1000	Mindset hardware model (16 channels)
MS24	Mindset hardware model (24 channels)
NCS	Nolan Computer Systems, L.L.C.
PC	Personal Computer
PCM	Pulse Code Modulation
RAM	Random Access Memory
RegEdit	Windows O/S Registry Edit program
O/S	Computer Operating System
SCSI	Small Computer Systems Interface
SPS	Samples per Second
SR	Sample Rate
V	Microvolts
V/pix	Microvolts per Pixel
VGA	Video Graphics Adapter card
X-axis	the horizontal axis
Y-axis	the vertical axis
Z-axis	the diagonal axis, representing <i>time</i>
\$	represents hexadecimal values

## Appendix B

### Using the Mindset Calibrator

#### B.1 INTRODUCTION

The *optional* Mindset24 Calibrator is a hardware device used in combination with the Mindset24 hardware and the MindMeld24 program. The calibrator injects a precise 16 Hz, 50 V signal simultaneously into all 24 Mindset channels. This calibration signal is read by the MindMeld24 LDC program and is used to calculate a correction factor (a *scaler*) for each channel to bring that channel to the signal norm.

Why is this necessary? Component manufacturing tolerances may cause amplifier gain to vary slightly among each channel and may cause the absolute gain of any individual amplifier to change with age. By injecting a known signal into the channels, a correction factor can be calculated for each channel to compensate for these changes. The scalers are saved in your PC's System Registry so that the MindMeld program can access and apply them.


If you intend to use Mindset with different PCs, access the Windows Registry Edit (RegEdit.exe) program to export the scaler data in a .reg file to your other PC's System Registry. The scaler data is in the *HKEY\_CURRENT\_USER* → *Software* → *NCS Programs* → *MindMeld24* folder. Consult with your local System/Network Administrator or Windows documentation for additional information on the RegEdit program. *Alternatively*, the Mindset calibration process may be performed with each PC.

It is recommended that you calibrate your Mindset once each year or whenever readings indicate that a calibration correction may be necessary.

*In this Appendix, a dollar sign (\$) is used to represent hexadecimal values.*

##### B.1.1 PARTS IDENTIFICATION

Your Mindset Calibrator is completely self-contained. The following illustration provides parts identification.



The Mindset Calibrator

##### B.1.2 THE CALIBRATION PROCESS

*It is assumed that at this point you have read the body of this documentation to familiarize yourself with the MindMeld applications, terminology and conventions of this document.*

###### B.1.2.1 Calibration Overview

A special calibration scheme named *calibrate.em3* is included on the floppy disk provided with the Mindset Calibrator. Use this scheme during calibration by copying it to your hard drive.

To calibrate your Mindset24, connect the calibrator to Mindset s expansion connector and run the special calibration scheme. Initially, this scheme provides verification that the calibrator is functioning properly. Secondly, it performs the actual calibration. When calibration is complete, quit the MindMeld24 application to save the calibration data into your PC s System Registry.

### B.1.2.2 The Calibration Process Instructions

#### CALIBRATION PROCESS INSTRUCTIONS


<b>User Action</b>	<b>Results/Comments</b>
<b><i>Warm-up Your Calibrator</i></b>	
<i>Turn on your Calibrator at least 20-to-30 minutes prior to conducting this calibration process.</i>	<i>The Calibrator must be completely warm before using. At this point in the process, the Calibrator can be turned-on without connecting it to your Mindset. (Or the Calibrator can be turned-on and connected to your Mindset now following the directions in the next two steps.)</i>
<b><i>Remove External Devices</i></b>	
<i>Remove any devices (Electro-cap, etc.) connected to your Mindset s input jacks.</i>	<i>When calibrating your Mindset, the Calibrator must be the only device connected during this process.</i>
<b><i>Connect Your Calibrator</i></b>	
<i>Connect your warmed-up Calibrator to Mindset s expansion connector. (If it is attached, first remove the Montage Selector.)</i>	<i>The Calibrator can be connected to your Mindset while power is on for both the Calibrator and for Mindset. (In other words, there is no need to turn off Mindset s power when connecting the Calibrator.)</i>
<b><i>The Calibration Process: Verify the Calibrator</i></b>	
<i>Start MindMeld24 LDC.</i>	<i>Refer to paragraph 2.3.</i>
<i>Load → Scheme → calibrate.em3.</i> <ul style="list-style-type: none"> <li>• If necessary, copy the special calibration scheme from the floppy disk to your hard drive.</li> <li>• If necessary, resize the MindMeld24 LDC main window and/or move the Stripchart tool window so you can view all 16 channels.</li> </ul>	<i>Refer to paragraph 4.1.1.</i>
<i>Start → Scheme.</i>	
<i>Observe the Stripchart tool.</i>	<i>Since an identical signal is injected into all 24 channels, the stripchart should appear as follows:</i>



- If the stripchart does not appear as in the above illustration, turn off the Calibrator's power for 5 seconds and turn it back on.
- If the stripchart still does not appear as in the above illustration, call NCS for technical assistance, or disconnect and return the Calibrator to your distributor or NCS.

Observe the Histogram tool. (If necessary, resize the MindMeld LDC main window and/or move the Histogram tool window to view the entire window.)

The histogram should appear as shown in the illustration below. (Refer to paragraph B.1.3 below for additional information.)



If your Stripchart and Histogram tool windows properly match these illustrations, *Stop* → *Scheme*.


At this point, your Calibrator appears to be functioning properly.

### *The Calibration Process: Calculate Calibration Scalars*

Select *Special* → *Calibrate*, the Calibration Procedures dialog box displays.

This *dialog* presents a summary of the entire calibration process. You have already completed the first three

items on this list.



If you are *satisfied* that the Calibrator is operating properly, click OK in the Calibration Procedures dialog box.

*or*

If you are *not satisfied* that the Calibrator is operating properly *or if you have not implemented any of the above steps*, click CANCEL in the Calibration Procedures dialog box.

If you clicked CANCEL, initiate or repeat as necessary the steps in this Calibration Process Instructions table.

If any problems persist, call NCS for technical assistance, or disconnect and return the Calibrator to your distributor or NCS.

*Start → Scheme.*

You are re-starting the special calibration scheme (calibrate.em3). The Calibration Window displays the calibration's progress. The meaning of the numbers in the Calibration Window are explained below. As the scheme runs, the word "Calculating..." shows in the far right column.

Calibration Window (averaged values, in hexadecimal)						
Channel	ADC Maximum	ADC Minimum	ADC Delta	Scaler	Scaled Delta	
1	98AD	46EA	5223	0.9CC	43F2	Calculating...
2	987E	443C	5746	0.847	43F3	Calculating...
3	97DE	46EC	5171	0.9CE	43F6	Calculating...
4	9849	4348	5800	0.84C	43F3	Calculating...
5	97FB	48E1	4F99	0.92E	43ED	Calculating...
6	98E1	48E7	5039	0.92C	43C8	Calculating...
7	9CB7	43EC	592A	0.82E	43D3	Calculating...
8	9D2F	44F2	583D	0.83E	43E8	Calculating...
9	9DE3	43E1	5A61	0.81C	4386	Calculating...
10	98FA	445F	5751	0.84F	4A1A	Calculating
11	9C37	4541	56F5	0.84F	4398	Calculating
12	994F	46FF	534F	0.85C	4A2F	Calculating
13	9754	46A6	508F	0.921	4A2F	Calculating
14	98C3	46F4	519F	0.9CF	4A18	Calculating
15	9807	43FF	5793	0.84F	4A76	Calculating
16	9576	475C	4C09	0.9F2	4A??	Calculating

Allow the scheme to run until "DONE!" shows next to all 24 channels.

Calibration Window (averaged values, in hexadecimal)						
Channel	ADC Maximum	ADC Minimum	ADC Delta	Scaler	Scaled Delta	
1	98A3	4671	5231	0.90C	43FE	DONE!
2	9B5C	43FA	575A	0.847	4A34	DONE!
3	97DE	46EE	5179	0.90E	43FC	DONE!
4	9BC4	43A0	5824	0.84C	4A11	DONE!
5	9815	48EE	4F49	0.92E	43FB	DONE!
6	987C	482E	5050	0.92C	43D3	DONE!
7	9CCD	43E7	5966	0.82E	4A39	DONE!
8	9D2B	44CD	585E	0.83E	4A34	DONE!
9	9DE4	43EF	5A54	0.82C	4A13	DONE!
10	9C2C	44D4	5758	0.84E	4A1F	DONFI
11	9C3F	454F	56FA	0.85E	4A78	DONFI
12	9973	464F	5333	0.89F	4A14	DONFI
13	977E	469F	508D	0.921	4A3E	DONFI
14	97F7	467C	51A5	0.9CF	4A1D	DONFI
15	9B1F	43FE	57A7	0.84E	4A1A	DONFI
16	9578	47F8	4F70	0.9E7	4A7A	DONFI

Stop → Scheme.	That's it! Your Mindset is calibrated. When you quit the MindMeld LDC application, the scalers are automatically saved into your PC's System Registry.
Quit and restart the MindMeld program to implement the calibration scalers.	

### B.1.2.3 Calibration Window Details

The columns in the Calibration Window (return to the figures in the Calibration Process Instruction table) are described below. Each horizontal row represents the data for one channel. A dollar sign (\$) is used to represent hexadecimal values.

#### Channel

The *hardware channel number* shows in this first column.

#### ADC Maximum and ADC Minimum

As the calibration signal is read from the Mindset hardware, these numbers are the *averaged peak (high and low) values* as reported by Mindset's 16-bit analog-to-digital converter (ADC). The no signal (zero) position is set to \$7000.

#### ADC Delta

This number is the *averaged difference between the ADC Maximum and ADC Minimum values*. Since there are rounding errors involved, the difference may not be exactly ADC Maximum minus ADC Minimum.

#### Scaler

For the Mindset hardware, a 50 V signal should produce an average (nominal) ADC Delta of \$4A25 plus or minus about \$95. If the calculated ADC Delta is above or below this nominal value, the scaler is adjusted to bring the ADC Delta to the nominal value.

## Scaled Delta

The Scaled ADC Delta is the ADC Delta value with the Scaler value applied. When the scaled data falls within the nominal range, the channel is calibrated and the word "DONE!" appears in the Progress column.


## Progress Column

This column indicates the status of the calibration:

- As calibration proceeds, the word "*Calculating*" appears next to each channel's information.
- When a channel has been calibrated, this text changes to "*DONE!*".
- If any channel is too far out of specification, the words "Out of range" appear to identify a faulty channel. If Out of range appears for any channel, your Mindset hardware should be returned to your distributor or NCS for repair.

### B.1.3 ANOMALOUS READINGS FROM THE CALIBRATOR

If the calibrator is not fully warmed up, you may observe some spectral energy in the frequencies around 16 Hz as shown in the following graphic:



Histogram Tool with Spectral Energy at 16 Hz


A small amount (5  $\mu$ V or less) of energy in frequencies around 16 Hz is to be expected. *If this energy exceeds about 10  $\mu$ V when the calibrator is fully warmed up, the calibrator requires adjustment and should be returned to your distributor or NCS.*

### B.1.4 APPLICATION OF SCALER DATA

*MindMeld LDC and AF use the scalers to correct the raw or saved EEG data prior to sending it to their tools. To demonstrate the value of using calibration, below are two graphics of a Hemispheric Tool window in MindMeld AF. Both windows display 50  $\mu$ V calibration data*


captured with a typical production Mindset MS-1000; the first graphic displays the data without scalers applied, the second graphic displays the data with scalers applied.

This first graphic (i.e., Hemispheric Tool Data *without* Scalers Applied) shows the data without using the scalers, *note the slight differences in channel gains. Special → Use Calibration Scalers* was not invoked.



Hemispheric Tool Data *without* Scalers Applied

This second graphic (i.e., Hemispheric Tool Data *with* Scalers Applied) shows the data using the scalers, *channel differences are corrected. Special → Use Calibration Scalers* was invoked.



Hemispheric Tool Data *with* Scalers Applied

## Appendix C

### Listing of Pre-configured Schemes

The following pre-configured schemes are included on your CD-ROM. You are encouraged to experiment with these schemes (you may wish to make a copy of the scheme first and modify the copy). Adjust tool parameters and observe the result(s). Add new tools to the scheme and connect them by clicking and dragging your mouse from tool to tool. See Section 5 for details.

These schemes assume Host ID 0, SCSI ID 2 for the Mindset24 hardware. If your configuration is different, open the Edit Scheme window, right-click on the Mindset24 tool to open its properties window, and select `Auto Scan`. This will automatically locate your Mindset24 hardware. Save the scheme to preserve the new settings. You can restore any scheme to its original settings by copying that scheme from the supplied CD-ROM to your hard drive.

`Histo&Neuro.em3` - Histogram & Neuromap tools. The streaming data are run through an FFT using a Blackman window.

- The Mindset tool is set to 256 samples/second with a SCSI block size of 768 bytes.
- The Neuromap tool is set with a threshold of 0 micro-volts and a range of 100 micro-volts.
- The Histogram tool is set to 24 channels, 2 Hz to 34 Hz and 50 micro-volt maximum displayed amplitude.

`Histogram.em3` - Histogram tool. The streaming data are run through an FFT using a Blackman window.

- The Histogram tool is set to 4 channels, 2 Hz to 34 Hz with x-axis labels on and a 50 micro-volt maximum displayed amplitude.

`Strip&CSA.em3` - Stripchart and Compressed Spectral Array (CSA) tools. Streaming data are sent to the CSA tool through an FFT using a Blackman window and averaged 8 times through the Averager tool.

- The Stripchart tool is set to 24 channels with an 8 second sweep and 5 micro-volts/pixel.
- The CSA tool is set to 19 channels, 2 Hz to 34 Hz and 50 micro-volt maximum displayed amplitude.

`Strip&Neuromap.em3` - Stripchart and Neuromap Tools. Streaming data are sent to the Neuromap tool through an FFT using a Blackman window.

- The Stripchart tool is set to 24 channels, 8 second sweep time and 5 micro-volts/pixel.
- The Neuromap tool is set to a threshold of 0 micro-volts and a range of 50 micro-volts.

`Stripchart.em3` - Stripchart Tool.

- The Stripchart Tool is set to display 24 channels with an 8 second sweep. The vertical sensitivity is set to 3 micro-volts/pixel.

`TenTwenty.em3` — Stripchart, Neuromap and Histogram Tools.

- The Stripchart Tool is set to display 19 channels with an 8 second sweep. The vertical sensitivity is set to 3 micro-volts/pixel. Each channel label is set to its standard 10-20 name.